

INVACARE® CANADA SCOOTERS



Yes, you can.®



Yes, you can.®

Invacare® scooters are designed by a team of engineers and clinicians, with input from casual and everyday riders. The result, is a family of scooters that check off every box in terms of performance, design and comfort.

Featuring everything from indoor micro-portable scooters, to outdoor robust, heavy duty scooters. The Invacare Canada scooter lineup is unparalleled in its offering and presents a solution for every need.

Tested at the world class TÜV certification facility in Germany, Invacare scooters have been certified for safety, function and quality for peace-of-mind reliability and long-lasting durability.






reddot award 2014
 winner

COLIBRI4B-12

Colibri™ 3-Wheel & 4-Wheel Micro-Portable

Model Numbers:

	3-Wheel 12 Ah Batteries	4-Wheel 12 Ah Batteries	4-Wheel 18 Ah Batteries
Orange	COLIBRI3O-12	COLIBRI4O-12	COLIBRI4O-18
Red	COLIBRI3R-12	COLIBRI4R-12	COLIBRI4R-18
White	COLIBRI3W-12	COLIBRI4W-12	COLIBRI4W-18
Green	COLIBRI3G-12	COLIBRI4G-12	COLIBRI4G-18
Blue	COLIBRI3B-12	COLIBRI4B-12	COLIBRI4B-18
Silver	COLIBRI3S-12	COLIBRI4S-12	COLIBRI4S-18

The Invacare® Colibri is the perfect travel companion. Designed with unparalleled durability features and enhanced with patented LiteLock™ Technology for simple disassembly, the Colibri is the ideal solution for trips to the mall, visits to the park or weekends away. Standard with batteries and a 300 lbs. weight capacity.

Features

- Three platforms with different wheel and battery options: 200 x 50 three-wheel with 12 Ah batteries, 200 x 50 four-wheel with 12 Ah batteries or 210 x 65 four-wheel with 18 Ah batteries.
- Patented LiteLock™ technology enables easy assembly and disassembly for transportation and travel.
- Six interchangeable colour panel options to personalize the Colibri to suit your unique style.
- Durable, double-spoke black wheels and composite front basket.
- Adjustable three-position tiller with delta steering handle.
- Padded vinyl driver seat featuring width and angle adjustable arms, height adjustment and stylish accent stitching.
- 12 Ah and 18 Ah dual battery pack options offer convenient charging through the front port and up to a 16 km travel range.
- Provider-programmable electronics to reverse throttle direction and adjust driving performance.
- Standard built-in safety features include rear anti-tippers, front and rear mounted reflectors and brake disengagement lever.
- 300 lbs. weight capacity.

STYLE IS FOREVER. INVACARE SCOOTERS.

Colibri™ Scooter Features



Six Colour Panel Options

Personalize your scooter with one of six interchangeable colour options available. Choose from Diamond Silver, Sunstone Orange, Jade Green, Jasper Red, Sapphire Blue or Marble White.



LiteLock™ Technology

Colibri's patented LiteLock system allows you to simply remove the battery box to disengage the locking mechanism for easy disassembly of the scooter into manageable components for travel.



Comfortable and Adaptable

The three-position tiller with delta steering handles ensure a comfortable driving position, while the padded drivers' seat features width and angle adjustable arms and stylish accent stitching.



Standard Safety Features

Rear anti-tippers, front and rear reflectors and a manual brake disengagement lever are standard on Colibri scooters. The reinforced frame enhances the durability of the scooter through daily travels.



Dual Options for Tires and Batteries

Sleek double-spoked black rims on non-marking low profile tires available in 200 x 50 and 210 x 65 sizes. 12 Ah and 18 Ah battery options provide convenient charging on and off the scooter and a travel range of up to 16 km.



Five Lightweight Components

- Seat 22 lbs
- Front 3W 29 lbs
- Front 4W 34 lbs
- Rear 21 lbs
- Battery Pack 12 Ah 22 lbs
- Battery Pack 18 Ah 32 lbs
- Basket 1.45 lbs

Product Specifications

Overall Dimensions	Colibri 3W 12Ah/ 4W 12Ah/4W 18Ah		Leo 3W/4W		Pegasus 3W/4W		Pegasus 4W		Comet 4W	
					METRO	PRO	PRO	ULTRA		
Base Width	19.25"	23.6"	26"	25.5"	26"	26"	26"	26"	26"	26"
Max. climbable obstacle height	1.75"/1.75"/3"	2.4"/2.5"	3"	4"	3"	4"	4"	4"	4"	4"
Length	40"	47.5"	49"/50"	52"	49"/50"	52"	57"	57"	57"	57"
Height	34.1"	36"	49.5"	49.5"	49.5"	49.5"	49.5"	50.75"	50.75"	50.75"
Usable Foot Space Length	18"/12"/12"	22"/13"	23"/13"	13"	23"/13"	13"	13.5"	13.5"	13.5"	13.5"
Seat Dimensions										
Seat-to-Floor	14.5" - 16.5"	16.5" - 19.5"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"
Width	17"	19"	20"	20"	20"	20"	20"	20"	20"	26"
Depth	17.5"	16"	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"	22.5"
Back Height	15.5"	18.5"	27" - 29"	27" - 29"	27" - 29"	27" - 29"	27" - 29"	27" - 29"	27" - 29"	32" - 35"
Back Angle	N/A	N/A	90° - 130°	90° - 130°	90° - 130°	90° - 130°	90° - 130°	90° - 130°	90° - 130°	95° - 130°
Arm Dimensions										
Length	11.5"	14"	14"	14"	14"	14"	14"	14"	14"	14"
Seat-to-Top of Arms	8"	8"	7" - 10"	7" - 10"	7" - 10"	7" - 10"	7" - 10"	7" - 10"	7" - 10"	5" - 7"
Width Between Arms	16" - 19"	23" - 27"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"
Recommended Battery										
	2 x 12 V, 12 Ah*/	2 x 12 V, 36 Ah*	2 x 12 V, 40 Ah	2 x 12 V, 75 Ah	2 x 12 V, 40 Ah	2 x 12 V, 75 Ah	2 x 12 V, 75 Ah	2 x 12 V, 75 Ah	2 x 12 V, 75 Ah	2 x 12 V, 75 Ah
	2 x 12 V, 18 Ah*									
Off-Board Charging System										
	2 Ah/2 Ah/5Ah	5 Ah	8 Ah	8 Ah	8 Ah	8 Ah	8 Ah	8 Ah	8 Ah	8 Ah

*included with scooter

Product Specifications

	Colibri 3W 12Ah/ 4W 12Ah/4W 18Ah	Leo 3W/4W	Pegasus 3W/4W METRO	Pegasus 4W PRO	Comet 4W PRO	Comet 4W ULTRA Pro
Speed	8 km/h (5 mph)	8 km/h (5 mph)	12 km/h (7.5 mph)	15 km/h (9.5 mph)	15 km/h (9.5 mph)	10 km/h (6 mph)
Estimated Travel Range*	11 km (7 miles) / 11 km (7 miles) / 16 km (10 miles)	36 km (22 miles)	42 km (26 miles)	55 km (34 miles)	55 km (34 miles)	58 km (36 miles)
Weight Limit (lbs.)	300**	350**	300**	350**	350**	485**
Maximum Incline Capability	6°	6°/8°	6°	10°	10°	10°
Tires						
Dimensions	200 x 50 / 200 x 50 / 210 x 65	9" x 3.5"	11" x 3.5" 4.10/3.50-5	12" x 3.5" 3.50-5	13" x 4" 4.00-5	13" x 4" 4.00-5
Tire Pressure	N/A	50 psi	50 psi	50 psi	50 psi	50 psi
Weight (lbs.)	91.6 / 97.4 / 108 (inc. batteries)	174/184 (inc. batteries)	114/125	130	138	155

Warranty Two-year limited warranty against manufacturing defects (excluding upholstered materials, padded materials and tires/wheels). Six months on batteries purchased from Invacare®. Manufacturer's warranty does not cover normal wear.

*Travel range is an estimate only based on ISO 7176-2 standards determined by testing using a 100 kg/220 lb. dummy. Travel range will vary depending on user weight, speed, battery type, battery condition and charge level, weather, terrain and condition of the drive wheels.

**Weight limit is the maximum total weight of the rider, accessories and additional items the rider may be transporting. Example: If the weight limit of the scooter is 350 lbs. and additional items equal 25 lbs., the maximum rider weight can be 325 lbs.



Yes, you can.[®]

Invacare Canada LP

www.invacare.ca

Canada

570 Matheson Blvd. E.

Unit 8

Mississauga, Ontario

L4Z 4G4 Canada

(800) 668-5324

All rights reserved.
Trademarks are identified
by the symbols TM, SM and [®].
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted.

Specifications are subject to
change without notification.

© 2018

Invacare Canada

Form No. 08-621C

Rev 06/18