

# INVACARE<sup>®</sup> CANADA SCOOTERS



Yes, you can.™



***Yes, you can.®***

Invacare® scooters are designed by a team of engineers and clinicians, with input from casual and everyday riders. The result, is a family of scooters that check off every box in terms of performance, design and comfort.

Featuring everything from indoor micro-portable scooters, to outdoor robust, heavy duty scooters. The Invacare Canada scooter lineup is unparalleled in its offering and presents a solution for every need.

Tested at the world class TÜV certification facility in Germany, Invacare scooters have been certified for safety, function and quality for peace-of-mind reliability and long-lasting durability.





## Leo™ 3-Wheel & 4-Wheel Indoor/Outdoor

Model Numbers:

	<b>3-Wheel</b>	<b>4-Wheel</b>
Blue	LEO-3B	LEO-4B
Red	LEO-3R	LEO-4R
Silver		LEO-4S

Available in three and four-wheel platforms, the Leo scooter is the ideal indoor/outdoor scooter due to its compact size and long-range capacity. Engineered with the highest regards to safety, function and quality, the Leo scooter is the optimal choice for those who value their independence. Features such as delta-steering handles, an integrated lighting system and a comfort-adjustable tiller and seat ensure that any journey on the Leo scooter is comfortable and safe. Includes 36 Ah batteries.

### Features

- Standard with 36 Ah batteries providing up to a 36 km range.
- Tool-free disassembly into five lightweight parts for travel and transportation.
- Compact size ensures easy control and maneuverability.
- Built-in splash guards protect the motor and electronics.
- Black, pneumatic tires are non-scuffing and provide a soft ride.
- 10" black pneumatic tires with polished aluminum rims.
- Two-step disengaging lever guards against unintentional disengagement out of drive mode.
- 360° rotating seat is height adjustable and features width and angle adjustable arms and stylish accent stitching.
- Full lighting system. Efficient LED lights increase visibility and improve scooter range.
- Provider-programmable electronics to reverse the throttle direction and adjust driving performance.
- Infinitely adjustable tiller with ergonomic delta-steering handles.
- Easy to use, lit control panel with programmable audible and visual signals.
- Dual rear accessory ports accommodate the use of multiple accessories at one time.
- 350 lbs. weight capacity.

## Leo™ Scooter Features

---



### Infinite Tiller Adjustment

The angle of the tiller easily adjusts to ensure a comfortable driving position and to allow for easy control panel interaction and viewing.



### Comfort Seating

Height adjustable seat features width and angle adjustable arms, fore/aft adjustment and 360° of rotation. Dual accessory ports allow for the use of two rear accessories, if required.



### Brake Light/Indicators

Rear brake lights engage upon the release of the throttle. Left and right direction indicators automatically stop blinking after 30 seconds.



### Intuitive Control Panel

Designed to ensure easy and accurate operation with domed buttons, LED lights and audible signals as well as a lit battery indicator showing remaining battery power.



### Splash Guards

Built-in splash guards protect the motor and electronics and help prevent unnecessary servicing.

## Product Specifications

Overall Dimensions	Colibri 3W 12Ah/ 4W 12Ah/4W 18Ah		Leo 3W/4W		Pegasus 3W/4W		Pegasus 4W		Comet 4W	
					METRO	PRO	PRO	ULTRA		
Base Width	19.25"	23.6"	26"	25.5"	26"	26"	26"	26"	26"	26"
Max. climbable obstacle height	1.75"/1.75"/3"	2.4"/2.5"	3"	4"	3"	4"	4"	4"	4"	4"
Length	40"	47.5"	49"/50"	52"	49"/50"	52"	57"	57"	57"	57"
Height	34.1"	36"	49.5"	49.5"	49.5"	49.5"	49.5"	50.75"	50.75"	50.75"
Usable Foot Space Length	18"/12"/12"	22"/13"	23"/13"	13"	23"/13"	13"	13.5"	13.5"	13.5"	13.5"
<b>Seat Dimensions</b>										
Seat-to-Floor	14.5" - 16.5"	16.5" - 19.5"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"
Width	17"	19"	20"	20"	20"	20"	20"	20"	20"	26"
Depth	17.5"	16"	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"	22.5"
Back Height	15.5"	18.5"	27" - 29"	27" - 29"	27" - 29"	27" - 29"	27" - 29"	27" - 29"	27" - 29"	32" - 35"
Back Angle	N/A	N/A	90° - 130°	90° - 130°	90° - 130°	90° - 130°	90° - 130°	90° - 130°	90° - 130°	95° - 130°
<b>Arm Dimensions</b>										
Length	11.5"	14"	14"	14"	14"	14"	14"	14"	14"	14"
Seat-to-Top of Arms	8"	8"	7" - 10"	7" - 10"	7" - 10"	7" - 10"	7" - 10"	7" - 10"	7" - 10"	5" - 7"
Width Between Arms	16" - 19"	23" - 27"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"
<b>Recommended Battery</b>										
	2 x 12 V, 12 Ah*/	2 x 12 V, 36 Ah*	2 x 12 V, 40 Ah	2 x 12 V, 75 Ah	2 x 12 V, 40 Ah	2 x 12 V, 75 Ah	2 x 12 V, 75 Ah	2 x 12 V, 75 Ah	2 x 12 V, 75 Ah	2 x 12 V, 75 Ah
	2 x 12 V, 18 Ah*			Group 24		Group 24	Group 24	Group 24	Group 24	Group 24
<b>Off-Board Charging System</b>										
	2 Ah/2 Ah/5Ah	5 Ah	8 Ah	8 Ah	8 Ah	8 Ah	8 Ah	8 Ah	8 Ah	8 Ah

\*included with scooter

## Product Specifications

	Colibri 3W 12Ah/ 4W 12Ah/4W 18Ah	Leo 3W/4W	Pegasus 3W/4W METRO	Pegasus 4W PRO	Comet 4W PRO	Comet 4W ULTRA Pro
<b>Speed</b>	8 km/h (5 mph)	8 km/h (5 mph)	12 km/h (7.5 mph)	15 km/h (9.5 mph)	15 km/h (9.5 mph)	10 km/h (6 mph)
<b>Estimated Travel Range*</b>	11 km (7 miles) / 11 km (7 miles) / 16 km (10 miles)	36 km (22 miles)	42 km (26 miles)	55 km (34 miles)	55 km (34 miles)	58 km (36 miles)
<b>Weight Limit (lbs.)</b>	300**	350**	300**	350**	350**	485**
<b>Maximum Incline Capability</b>	6°	6°/8°	6°	10°	10°	10°
<b>Tires</b>						
<b>Dimensions</b>	200 x 50 / 200 x 50 / 210 x 65	9" x 3.5"	11" x 3.5" 4.10/3.50-5	12" x 3.5" 3.50-5	13" x 4" 4.00-5	13" x 4" 4.00-5
<b>Tire Pressure</b>	N/A	50 psi	50 psi	50 psi	50 psi	50 psi
<b>Weight (lbs.)</b>	91.6 / 97.4 / 108 (inc. batteries)	174/184 (inc. batteries)	114/125	130	138	155

**Warranty** Two-year limited warranty against manufacturing defects (excluding upholstered materials, padded materials and tires/wheels). Six months on batteries purchased from Invacare®. Manufacturer's warranty does not cover normal wear.

\*Travel range is an estimate only based on ISO 7176-2 standards determined by testing using a 100 kg/220 lb. dummy. Travel range will vary depending on user weight, speed, battery type, battery condition and charge level, weather, terrain and condition of the drive wheels.

\*\*Weight limit is the maximum total weight of the rider, accessories and additional items the rider may be transporting. Example: If the weight limit of the scooter is 350 lbs. and additional items equal 25 lbs., the maximum rider weight can be 325 lbs.



*Yes, you can.<sup>®</sup>*

**Invacare Canada LP**

[www.invacare.ca](http://www.invacare.ca)

**Canada**

570 Matheson Blvd. E.

Unit 8

Mississauga, Ontario

L4Z 4G4 Canada

(800) 668-5324

All rights reserved.

Trademarks are identified  
by the symbols <sup>TM</sup>, <sup>SM</sup> and <sup>®</sup>.

All trademarks are owned  
by or licensed to Invacare  
Corporation unless  
otherwise noted.

Specifications are subject to  
change without notification.

© 2018

Invacare Canada

Form No. 08-621C

Rev 06/18