

SPECIFICATIONS


Product Weight:	12 kg (26.5 lbs.) ¹
Weight Capacity:	113 kg (250 lbs.) (159 kg (350 lbs.) with HD configuration)
Seat-to-Floor Heights:	38 cm - 49.5 cm (15"-19.5") front 36.8 cm - 48.3 cm (14.5"-19") rear 31.8 cm - 42 cm (12.5"-16.5") front (Ultra Hemi) 31.8 cm - 48.3 cm (12.5"-19") rear (Ultra Hemi)
Seat Sizes:	Width: 35.5 cm - 53.3 cm (14"-21") (HD: 45.7 cm - 56 cm (18" - 22")) Depth: 35.5 cm - 50.8 cm (14"-20") (HD: 45.7 cm - 50.8 cm (18" - 20")) (Available in 2.5 cm (1") increments)

COLOURS

METALLIC



Cosmic Silver Pacific Blue Plum Crazy



Purple Haze Rainforest Green

NON-METALLIC



Yellow Sunshine Blue Sky Green Mantis Pink Flamingo



Red Hot Black Abyss Orange Crush Arctic White

OPTIONS

REAR WHEELS:	51 cm, 56 cm, 61 cm, 63.5 cm & 66 cm (20", 22", 24", 25" & 26") Spinergy wire spoke 56 cm, 61 cm, 63.5 cm & 66 cm (22", 24" 25" & 26") Spinergy Light Extreme (available in red, yellow, blue, white, pink, green, orange and black) 61 cm, 63.5 cm & 66 cm (24", 25" & 26") Spinergy Sport Light Extreme (available in red, yellow, blue, white, pink, green, orange and black) 51 cm, 56 cm & 61 cm (20", 22" & 24") Composite mag 61 cm (24") heavy duty Skyway Tuff composite mag wheels Magic Wheels One-arm drive system		
TYRES:	56 cm, 61 cm, 63.5 cm & 66 cm (22," 24", 25" & 26") 110 PSI Primo tyre 61 cm & 66 cm (24" & 26") Aggressive knobby tyre 51 cm, 56 cm, 61 cm & 66 cm (20", 22", 24" & 26") Everyday grey tyre 51 cm, 56 cm, 61 cm & 63.5 cm (20", 22", 24" & 25") Polyurethane tyre (61 cm (24") also available in light grey) 66 cm (26") Kenda mountain bike knobby tyres (available in black or grey) 61 cm & 66 cm (24" & 26") 110 PSI Primo coloured tyres (available in red, blue or yellow) 56 cm, 61 cm & 63.5 cm (22", 24", 25" & 26") Schwalbe Marathon Plus and Marathon Plus Evolution		
HANDRIMS:	51 cm, 56 cm, 61 cm, 63.5 cm & 66 cm (20", 22", 24", 25" & 26") Aluminum 56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") Vinyl-coated 51 cm, 56 cm, 61 cm, 63.5 cm & 66 cm (20", 22", 24", 25" and 26") Projection (vertical or 45° oblique) 56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") Natural-Fit (available in standard or super grip) 61 cm (24") Natural-Fit LT (available in standard or super grip) 61 cm (24") Spinergy FlexRims 51 cm, 56 cm, 61 cm (20", 22" & 24") Chrome		
AXLE PLATES:	Standard (0°, 1.5°, 3° camber options) Amputee (0°, 1.5°, 3° camber options)		
CASTERS:	10 cm, 13 cm, 15 cm & 20 cm (4", 5", 6" & 8") composite 10 cm, 13 cm & 15 cm (4", 5" & 6") soft roll 13 cm (5") aluminum 10 cm (4") light-up casters 10 cm (4") coloured casters (available in red, blue, green, orange, yellow or clear) 10 cm, 13 cm, 15 cm & 20 cm (4", 5", 6" & 8") black, smooth 3-spoke tyres 15 cm (6") pneumatic and 20 cm (8") pneumatic or flat-free polyurethane casters		
FORKS:	13 cm (5") or 20 cm (8") quick release forks 13 cm (5") or 20 cm (8") bolt-on forks Frog Legs suspension forks 60 mm forks (for use with 20 cm (8") pneumatic casters only) Caster pin locks		
WHEEL LOCKS:	Pull-to-lock Push-to-lock Scissor locks	Pull-to-lock with extensions Push-to-lock with extensions Hill holder	Cable-linked, unilateral wheel locks
ARMRESTS:	30.5 cm (12") tubular swing-away (desk length) Height-adjustable drop-in (desk or full length) Locking or non-locking, flip-up (depth adjustable only) Height-adjustable, flip-back, dual post armrests (desk or full length)		
BACK CANES:	Straight, with or without push handles Angle-adjustable, with or without push handles 8 degree with push handles Straight reinforced Depth-adjustable frame option with fold forward back canes		
BACKS:	Standard nylon sling with convenience pocket Tension-adjustable upholstery Synergy® Shape Back		
FRONT RIGGINGS:	Aluminum, bi-directional swing-away leg rests, 60°, 70°, 70° tapered, 80° and 90° Elevating leg rests with calf pads Comfort adjust footplates Angle-adjustable foot plates Nylon heel loops		
ACCESSORIES:	Height-adjustable, removable anti-tip wheels Composite side guards Aluminum side guards Lap belt Quantum® backpack Spoke protectors Transit securement package Gel tubular arm cover Case Logic underseat organizer bag Combat, cane and crutch holder Case Logic cup holder		

¹) Standard configuration: 40.6 cm x 40.6 cm (16" x 16"), 48.3 cm (19") backrest, sling upholstery, push-to-lock brakes, 61 cm (24") spoke wheels, 3.5 cm (1-3/8") pneumatic tyres, 15.2 cm (6") casters (not including leg rests)

NOTE: Not all options available for all configurations. Please check order form for compatibility.

Litestream[®] XF



“Experience the performance and durability that I demand.”

– Bryan Anderson
Quantum Rehab® Spokesperson



20-24 Apollo Drive, Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

Tel (03) 8770 9600 • Fax (03) 9703 2720 • www.pridemobility.com.au

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

QUANTUM
LIFE BEYOND LIMITS™

www.quantumrehab.com.au



The Litestream® XF

Ultralight manual wheelchair uses a unique, innovative, rigidising double X folding frame that delivers the versatility and performance of a rigid frame while folding to a width of only 23 cm (9") for easy transport. The Litestream XF accepts a broad selection of accessories to meet a wide variety of needs and personal preferences while remaining an excellent overall value. The Litestream XF is also available in an HD configuration that weighs just 18 kg (39 lbs.) and features a cane-to-cane rigidising bar for added stability.

REAR WHEELS

CASTERS

TYRES



51 cm, 56 cm, 61 cm, 63.5 cm & 66 cm (20", 22", 24", 25" & 26") high-performance Spinergy wheels feature double-wall rims with straight-pull spokes for reduced weight and optimal durability

56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") Spinergy Light Extreme 12 spoke PBO Fiber wheels

AVAILABLE COLOURS:



61 cm, 63.5 cm & 66 cm (24", 25" & 26") Spinergy Sport Light Extreme 24 spoke PBO Fiber wheels

AVAILABLE COLOURS:



51 cm, 56 cm & 61 cm (20", 22" & 24") Skyway Tuff composite mag wheels are exceptionally durable and low maintenance

Also available:

61 cm (24") heavy duty Skyway Tuff composite mag wheels with aluminum handrim



Magic Wheels feature two gears to address real-life situations – 1st gear, for level ground and 2nd gear, providing a 100-percent increase in hill climbing torque



61 cm (24") lightweight, quick-release, one-arm drive system with aluminum handrim



10 cm, 13 cm, 15 cm & 20 cm (4", 5", 6" & 8") composite casters feature precision bearings and long-life poly tyres



10 cm, 13 cm & 15 cm (4", 5" & 6") soft-roll casters smooth lumps and bumps in everyday use



13 cm (5") aluminum casters with precision bearings and polyurethane tyres



10 cm (4") light-up casters



10 cm (4") coloured casters available in red, blue, green, yellow, orange or clear



10 cm, 13 cm, 15 cm & 20 cm (4", 5", 6" & 8") black, smooth 3-spoke tyres



15 cm (6") pneumatic casters



20 cm (8") pneumatic or flat-free polyurethane casters

HANDRIMS



51 cm, 56 cm, 61 cm, 63.5 cm & 66 cm (20", 22", 24", 25" & 26") Projection Handrim (vertical or 45° oblique)



61 cm (24") Spinergy FlexRims



56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" and 26") Natural-Fit LT Handrim (available in standard or super grip)

(Aluminum, vinyl-coated & chrome handrims also available but not shown)



61 m & 66 cm (24" & 26") Aggressive Knobby Tyre



51 cm, 56 cm, 61 cm & 66 cm (20", 22", 24" & 26") Everyday Grey Tyre



51 cm, 56 cm, 61 cm & 63.5 cm (20", 22", 24" & 25") Polyurethane Tyre (61 cm (24") also available in light grey)




56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") 110 PSI Primo Tyre



61 cm & 66 cm (24" & 26") 110 PSI Primo Coloured Tyre



66 cm (26") Kenda Mountain Bike Knobby Tyres



56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") Schwalbe Marathon Plus Evolution



56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") Schwalbe Marathon Plus

FORKS



Frog Legs suspension fork



20 cm (8") fork



13 cm (5") fork



Quick-release 60 mm fork (for use with 20 cm (8") pneumatic or flat-free polyurethane casters only)



Caster pin lock (available as bolt-on only)

WHEEL LOCKS



Pull-to-lock



Scissor lock



Pull-to-lock with extension



Push-to-lock with extension



Hill holder



Cable-linked, unilateral wheel locks (for use with one-arm drive system)



Push-to-lock

LEG RESTS



Elevating bi-directional leg rests with swing-away calf pads



60°, 70°, 70° tapered, 80° and 90° aluminum bi-directional swing-away leg rests



10-degree comfort adjust footplates in 11 cm, 13 cm & 15 cm (4.25", 5" & 6") widths




Angle-adjustable footplates in 14 cm, 16.5 cm & 19 cm (5.5", 6.5" & 7.5") widths



Nylon heel loops

BACK OPTIONS



Depth-adjustable frame option with fold forward back canes



Tension-adjustable upholstery



Standard nylon sling



Synergy® SHAPE back

TRANSIT SECUREMENT



The Transit Securement Package allows convenient securement of the Litestream® XF for safe transport.



ACCESSORIES AND PERSONALISATION



Height-adjustable anti-tips



Composite side guards



Aluminum side guards



Spoke protector



Case Logic cup holder



Quantum® backpack



Case Logic underseat organizer bag



Combat, cane & crutch holder

ARMRESTS



30.5 cm (12") tubular, swing-away, desk-length armrests (shown with optional gel arm covers)



Height-adjustable, flip-back, dual post armrests (desk or full length)



Locking or non-locking, flip-up armrests (for use with depth-adjustable option only)



Height-adjustable drop-in armrests (desk or full length)