

VivaLIFT![®]
POWER RECLINERS
By Pride[®] Mobility Products Corp.

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER

WITH
COUNTLESS
POSITION
COMBINATIONS



Shown in Stonewash
Granite Fabric

URBANA



Innovative technology meets contemporary style for a recliner that supports you and your busy life. The Urbana Collection by **VivaLift!**[®] Power Recliners is the perfect addition to your townhouse or condo. After a long day, experience the comfort of power lumbar, which can be adjusted at the touch of a button. The Urbana lets you fully recline with infinite, lay-flat positioning and experience maximum relaxation.

Pride Mobility[®] (US) 800-800-8586 | (CAN) 888-570-1113
pridemobility.com | prideliftchairs.com |    

URBANA

COLLECTION	URBANA v.1
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift
WEIGHT CAPACITY	375 lbs.
SEAT-TO-FLOOR	20"
SEAT DEPTH	20"
SEAT WIDTH	21"
TOP OF BACK TO SEAT	25"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	42"
DISTANCE FROM WALL	14"
OVERALL WIDTH	33"
CHAIR WEIGHT	145 lbs.
WARRANTY	7-Yr. Pro-rated


FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

POSITION GUIDE:

-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  INFINITE LAY FLAT

STANDARD FABRIC: STONEWASH



GRANITE

GUNMETAL

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.



Pride Mobility® (US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    