COLLECTION

## **Oasis** COLLECTION



















**POWER LIFT RECLINER** 

**END-TO-END** 



	-				A CONTRACTOR	1	1		1
MODEL	LC-380	LC-580iM	LC-580iL	LC-525iPW	LC-525iS	LC-525iM	LC-525iL	LC-108	LC-515iL
Size	Medium/Large	Medium	Large	Petite Wide	Small	Medium	Large	Medium	Large
Position	3-position Chaise Lounger 5'4" - 6'	True Infinite-position, Trendelenburg Chaise Lounger 5'4" - 6'	True Infinite-position, Trendelenburg Chaise Lounger 5'9" - 6'2"	True Infinite-position, Trendelenburg Chaise Lounger 5'2" and below	True Infinite-position, Trendelenburg Chaise Lounger 5'3" and below	True Infinite-position, Trendelenburg Chaise Lounger 5'4" - 6'	True Infinite-position, Trendelenburg Chaise Lounger 5'9" - 6'2"	Infinite-position, zero gravity, Trendelenburg Chaise Lounger 5'4" - 6'	True Infinite-position, Trendelenburg Zero Gravity 5'4" - 6'
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	21"	20"	21"	19"	19"	20"	21"	19"	20"
Seat Depth	21"	20"	21"	19"	19"	20"	22"	19"	20"
Seat Width	24.5"	20"	21"	24"	19"	20"	22"	20"	21"
Top of Back to Seat	25"	26"	29"	28"	26"	29"	30"	29"	25"
Back Style	Blown fiber biscuit back	Blown fiber biscuit back (lumbar pillow standard)	Blown fiber biscuit back (lumbar pillow standard)	Blown fiber split back (lumbar pillow standard)	Blown fiber split back (lumbar pillow standard)	Blown fiber split back (lumbar pillow standard)	Blown fiber split back (lumbar pillow standard)	Sewn pillow back (back upholstery is fully removable and adjustable) (lumbar pillow standard)	Blown fiber, Pillow back
Fully Padded Chaise	YES	YES	YES	YES	YES	YES	YES	YES (seat pad is fully removable and adjustable)	YES
Heat and Massage <sup>3</sup>	YES	YES	YES	YES	YES	YES	YES	NO	YES
Additional Options	Footrest extension (std), heat and massage, wireless remote	Deluxe heat and massage	Deluxe heat and massage	Deluxe heat and massage	Deluxe heat and massage	Deluxe heat and massage	Deluxe heat and massage	Mattress option, footrest extension (standard)	Footrest extension and lumbar pillow (standard)
Seat Type	Pocketed coil springs w/ memory foam topper	Pocket coil seat with memo- ry foam topper	Pocket coil seat with memory foam topper	Pocketed coil springs	Pocketed coil springs	Pocketed coil springs	Pocketed coil springs	Sinuous spring	Pocketed coil springs
Distance from Wall	20"	20"	21"	19"	20"	20"	21"	24"	21.5"
Overall Width	40"	39"	39"	38"	34"	34.5"	37"	31.5"	33.5"
Chair Weight	120 lbs.	130 lbs.	132 lbs.	130 lbs.	130 lbs.	130 lbs.	130 lbs.	136 lbs.	131 lbs.
Hand Control	Easy Touch with USB Port	Easy Touch Infinite- Position with USB port	Easy touch Infinite-position with USB port	Easy touch Infinite- position with USB port	Easy touch Infinite- position with USB port	Easy touch Infinite- position with USB port	Easy touch Infinite- position with USB port	Easy touch Infinite- position with USB port	Easy touch Infinite-position with USB port
Warranty	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated
Standard  Optional			TIMBER OAT DEEP EMBER				Cloud 9  (Quick Ship applies to LC-515L only) (Spruce not available on LC-515L)		
Optional	Crypton Aria®	Ultraleather <sup>™</sup>		Crypton Aria®	Ultraleather"	1	Lexis Sta-Kleen®	Lexis Sta-Kleen°	Ultraleather <sup>™</sup>

BLACK PACIFIC SPRUCE STONE WALNUT







Pride





- 1. Unique spring design and a one-piece soft foam seat for outstanding relaxation.
- 2. Features a dual motor sync system.
  3. USB hand control not available with heat and massage options.

BLACK PACIFIC SPRUCE STONE WALNUT

- 4. Not all hand controls are available for all models. Contact customer service for applicable options.
- 5. Toggle hand control only available for single motor chairs.

Like all furniture, Pride® Power Lift Recliners are individually built and there may be differences in actual measurements due to variations in fabric, upholstery, and filling, in addition to the method of

NORTH AMERICA'S #1

POWER LIFT RECLINER




## PARTNER TRADEMARKS









REMOTE

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice

Beamsville, Ontario • 888-570-1113





COLLECTION

HOMEDÉCOR
















































LC-570M LC-570L I	LC-250 LC-102 LC-106	LC-220	LC-107
Medium Large	Medium Super Petite Medium	Duty Medium	Medium
sistion 3-position Infinite Chaise Lounger Chaise Lounger 5'9" - 6'2" 5'9" - 6'2"	3-position Chaise Lounger 5'4"- 6'0"  3-position Chaise Chaise Chaise Lounger Lounger 4'8" - 5'4"  5'2"- 5'10"	- 2-position C	True Infinite-position, Trendelenburg Chaise Lounger 5'4"- 6'0"
lbs. 375 lbs. 325 lb	375 lbs. 250 lbs. 325 lbs.	375 lbs. 3	325 lbs.
21" 19"	20" 16" 18"	20" 20	20"
22" 20"	20" 18" 19"	20"	20"
22" 19"	20" 16" 20"	20" 20	20"
26"	26" 23.5" 27.5"	24.5" 26	27.5"
vn fiber Blown fiber erfall back waterfall back Blown back	Blown fiber split back Dual pillow Sewn pillow back		Sewn pillow back (lumbar pillow standard)
YES YES	YES YES YES	e, fully NO Y	YES
YES NO	NO NO NO	NO N	NO
trest First Up™, footrest extension (standard), heat and sage footbase extension (standard), heat and (standard) (standa	N/A N/A	N/A N	Footrest extension (standard)
l Comfort <sup>1</sup> Total Comfort <sup>1</sup> Total C	Total Comfort <sup>1</sup> Sinuous spring Sinuous spring	Sinuous spring To	Sinuous spring
24" 24"	20" 16" 18"	11" 20	26"
38" 33"	34" 25.5" 32.75"	32.5" 34	35"
lbs. 114 lbs. 147 lb	118 lbs. 96 lbs. 96 lbs.	114 lbs. 11	116 lbs.
y touch with Easy touch with USB port USB p	Easy touch Easy touch	USB Easy touch E	Easy touch Infinite-position with USB port
Pro-rated 7-Yr. Pro-rated 7-Yr. F	2-Yr. Limited 1-Yr. Limited 1-Yr. Limited	2-Yr. Limited 2-	1-Yr. Limited
s to LC-450C only) vailable fabrics for LC-570P.	Cloud 9  BLACK PACIFIC SPRUCE STONE WAR	Cloud 9  BLACK PACIFIC SPRUCE CHERRY	ALNUT
ESPRE Dura		m r Lift Recliners.)	
	TIMBER OAT DEEP SKY  POSITION GUIDE (The lifted, seated and reading positions are achievable on all Pride® Power Lift Recliners.)	POSITION GUIDE (The lifted, seated and reading positions are achievable on all Pride® Powe	POSITION GUIDE (The lifted, seated and reading positions are achievable on all Pride® Power Lift Recliners.)



NORTH AMERICA'S #1 POWER LIFT RECLINER

## ADDITIONAL NOTES:

- All chairs have head and arm covers standard (Coming soon for Home Décor Collection)
- All chairs have an engineered furniture-grade laminate frame
- · Lifetime Limited on lift, recline and scissor mechanism and wood frame on all power lift recliners (Exceptions: Essential Collection (1-year), LC-805 (2-year), LC-358XL and LC-358XXL (7-year))
- Transformer features self-diagnostic electronics with 9' cord
- All chairs have integrated emergency battery backup (batteries not included)
- · All chairs can have fiber fill added or removed for custom comfort adjustment
- All chairs are shipped in an impact-resistant, reinforced, corrugated box



LIFTED





















- The 3-position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing.
- The infinite-position lay flat recline is perfect for taking a nap. Infinite-position Pride® Power Lift Recliners enable the backrest and footrest to move independently from one another so that a variety of positions can be reached. The LC-107, LC-580i and LC-525i can have the footrest fully elevated while causing minimal impact to the back angle, achieving true infinite adjustability. True infinite positioning allows for easier social engagement, eating and crafting activities while having the footrest fully elevated. Infinite-Position

The 2-position TV recline is achievable on all Pride® Power Lift Recliners. This is a comfortable position for performing many relaxing activities including reading, watching TV, working on your tablet or for conversation with family and friends.

The Trendelenburg position, elevates one's feet higher than his or her head for increased comfort and relaxation.

The zero gravity position, which can be reached on Infinite-position Pride® Power Lift Recliner models, positions one's thighs and lower legs higher than his or her torso.

LIVE YOUR BEST pridemobility.com